

CLINELL UNIVERSAL WIPES

- HOSPITAL GRADE DISINFECTING CLEANING WIPES



Clean and disinfect surfaces regularly with Clinell Universal Wipes. With a textured weave these wipes not only disinfect surfaces but offer a cleaning function too. Follow the instructions outlined below for the most effective way to clean and disinfect surfaces.

Clinell Universal Wipes are proven effective against MERS-CoV, a closely-related strain of virulent human coronavirus, according to EN14476.

1

Wipe in an "S" shaped pattern from top to bottom and clean to dirty to cover the entire area.

2

Once the wipe is dry or soiled, throw it out and replace with a new one.



3

Allow to air dry. Clinell Universal kills human coronavirus within 60 seconds.

4

Do not rewet the wipe or use the same wipe on multiple pieces of equipment.

HOW TO ORDER:

Stock Code: CW200

Pack Size: 200 wipes

MOQ: 1 pkt

Price: \$38.00+gst

DO NOT FLUSH DOWN TOILETS! DO NOT DISPOSE OF IN MACERATORS!